

Palintest

Water Analysis Technologies

Phosphate/12P

Photometric test for Total Phosphate in natural water and wastewater

COLOUR CHANGE

Colourless to Blue

RANGE:

0 - 12 mg/L P

TECHNICAL RESOURCE:

Phosphate 12P Test Method, Technical Information

01




Add 2 mL of sample to a Tubetest.

02




Use this to "blank" the instrument. Fit the Tubetest adapter to Lumiso, insert tube and press Blank.

03



Continue with the same tube by adding a Tubetest Phos No1 tablet. Crush and stir, until fully dissolved.

04



Then, add a Tubetest Phos No2 tablet and crush and stir.

05



Ensure cap is tight, then invert tube a few times to mix.

06



Place in the cell holder and press Measure.

07



Allow timer to count down 10 minutes. Concentration of Phosphate will then be displayed.