

Palintest

Water Analysis Technologies

Silica HR

Photometric test for Silica (high range) in natural and treated water

COLOUR CHANGE

Colourless to Yellow

RANGE:

0 - 150 mg/L SiO_2

TECHNICAL RESOURCE:

Silica HR Test Method, Technical Information

01




Fill to the 10 mL line with sample.

02




Use this to "blank" the instrument. Place in cell holder and press Blank.

03



Prepare sample by adding a Silica HR No1 tablet and crush and stir.

04



Then, add a Silica HR No 2 tablet and crush and stir.

05



Place in the cell holder and press Measure.

06



Allow timer to count down 10 minutes. The concentration of Silica will then be displayed.